

# AEROSEPT 500M

Airborne  
disinfection  
device



- Treats a room up to 300m<sup>3</sup>
- Complies with NF T 72-281 (2014 version)
- Fully automatic process with touch operator panel and delayed start
- For floors and surfaces disinfection with suitable associated products
- User friendly

## INDICATIONS

Standalone automatic atomiser, for airborne disinfection of previously cleaned surfaces.

To use with ASEPTANIOS AD (active from 30 minutes).

## ENHANCED PERFORMANCE PROCESSES

- Solution heated to a constant temperature of 55°C throughout the duration of spraying.
- Concentration perfectly stabilized at 7ml/m<sup>3</sup>

\* Refer to the Scientific File to consult all claims and studies.

# A EROSEPT 500M

## Airborne disinfection device

### SIMPLE PROCESS

- Fully automated process
- Delayed start preset to 1 min
- Delayed start with
- Easy replacement of the container without risk of splashing

### SECURE PROCESS

- **Stops automatically at the end of the cycle:** displays an indication that the cycle was successful with a request for validation by the operator.
- **No splashing:** container equipped with a tube to avoid direct contact with the product.

### INSTRUCTIONS FOR USE

- Place the unit in a corner of the room
- Automatic calculation of diffusion time
- Indication of the time the cycle will end
- Delayed start allowing the operator to leave the room before the treatment starts
- Automatic stop at the end of the programme
- Ventilate the room



On



Non-compliant cycle



Compliant cycle

### TECHNICAL CHARACTERISTICS

- Storage capacity: 5L container
- Empty weight: ≈ 30 Kg
- Nominal flow-rate: 0,9L/hour
- Dimensions: 900 H x 350 L x 400 W
- Power supply: 100/240 V – 50/60 Hz – 4?1 A
- Multi-languages: French, English, Spanish, Polish, Turkish, Portuguese, Italian

### PRECAUTIONS FOR USE

Before using the device, read the instructions. The products used in this device are Dangerous - Observe the precautions on the label of the product used (Drawn up according to European rules on classification and labelling of chemicals). Use biocides safely. Please read the label and product information before use.

### REFERENCES

AEROSEPT 500M..... 418.100  
ASEPTANIOS AD (4 cans de 5L)..... 2122.795

**Healthline Pharmaecticals Pvt Ltd.**  
**Office No. 402 Al Hafeez Heights**  
**Gulberg 3, Lahore, Pakistan**  
**Contact: +92 321 4085989, +92 42 36400774**  
**Email: info@healthlineint.com**